




Oxford BTD

Bottle Top Dispenser



Oxford Lab Products
Bottle Top Dispenser
Product Overview



Oxford Lab Products

Bottle Top Dispenser

Product Overview

Build for Precision

Made of highly chemical-resistant materials, such as PTFE and PFA, the Oxford Lab Products Bottle-Top Dispenser is the ideal choice for dispensing bases, acids and solvents.

Eco-friendly - while discharging the tubes in the re-circulation mode, in fast priming, the reagent is being redirected back into the bottle, eliminating any waste of precious reagents, ensuring safety and savings.

Unique piston mechanism - cleaning of piston and cylinder without affecting the calibration.

Dual graduation marks - readability of the volume from both sides
Springless valve design - leads to smooth functioning and high chemical resistance.

Easy handling - discharge tube with 360° rotation, makes the bottle label visible all the time / telescoping filling tube compatible with different sizes of reagent bottles / ergonomic grip.

Autoclavable at 120 degrees Celsius for 20 minutes.

Easy to calibrate - complies with ISO 8655-5 standard.

Adapters- comes with 4 bottle adapters for the most common bottle sizes.



Bottle Adapter Sizes



45 mm



40 mm



38 mm



28 mm

Oxford Lab Products
Bottle Top Dispenser
Ordering Information

Oxford Bottle Top Dispenser

Model#	Description	Volume
AD-2.5	AccuTop Bottle Top Dispenser, Without Re-circulaton	0.25-2.5ml
AD-05	AccuTop Bottle Top Dispenser, Without Re-circulaton	0.5-5ml
AD-10	AccuTop Bottle Top Dispenser, Without Re-circulaton	1-10 ml
AD-30	AccuTop Bottle Top Dispenser, Without Re-circulaton	2.5-30 ml
AD-60	AccuTop Bottle Top Dispenser, Without Re-circulaton	5-60 ml
ADR-2.5	AccuTop Bottle Top Dispenser, WithRe-circulaton	0.25-2.5 ml
ADR-05	AccuTop Bottle Top Dispenser, WithRe-circulaton	0.5-5ml
ADR-10	AccuTop Bottle Top Dispenser, WithRe-circulaton	1-10 ml
ADR-30	AccuTop Bottle Top Dispenser, WithRe-circulaton	2.5-30 ml
ADR-60	AccuTop Bottle Top Dispenser, WithRe-circulaton	5-60 ml